

## Sample Collection Tubes

### I. APPLICATIONS

1. Sample collection
2. Sample transportation
3. Sample storage
4. Dispensing samples

### II. TUBE FEATURES

1. Material: Strong, medical-grade polypropylene, high-density polyethylene

### III. TUBE DESIGN

1. After sample collected, put the cotton swab into the tube, and tighten the cap, then the cotton swab will be clipped into the cap automatically.
2. Caps are external threads which will promote secure sealing.
3. 6 colors for caps to identify samples easily.
4. Inert PP and PE materials are resistant chemical corrosion.
5. Tubes are clear for convenient observation of samples and tests.

### IV. PACKAGING

Item	Description	Packaging	Sterile
10-0052	Sample collection tubes without caps, 6ml, PP, bulk	4000 Tubes/Case	No
10-0000	Caps for 10-0052 tubes, PE, white	4000 Caps/Case	No
10-0001	Caps for 10-0052 tubes, PE, red	4000 Caps /Case	No
10-0002	Caps for 10-0052 tubes, PE, green	4000 Caps /Case	No
10-0003	Caps for 10-0052 tubes, PE, blue	4000 Caps /Case	No
10-0004	Caps for 10-0052 tubes, PE, yellow	4000 Caps /Case	No
10-0005	Caps for 10-0052 tubes, PE, violet	4000 Caps /Case	No
10-0082	Sample collection tubes without caps, 12ml, PP, bulk	2100 Tubes/Case	No
10-0100	Caps for 10-0082 tubes, PE, white	10500 Caps /Case	No
10-0101	Caps for 10-0082 tubes, PE, red	10500 Caps /Case	No
10-0102	Caps for 10-0082 tubes, PE, green	10500 Caps /Case	No
10-0103	Caps for 10-0082 tubes, PE, blue	10500 Caps /Case	No
10-0104	Caps for 10-0082 tubes, PE, yellow	10500 Caps /Case	No

10-0105	Caps for 10-0082 tubes, PE, violet	10500 Caps /Case	No
---------	------------------------------------	------------------	----

## V. STERILITY


Not sterile and sterile tubes can also be ordered. For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

## VI. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


## VII. PRODUCT PARAMETER

Dimensions=mm



10-000X

10-010X



10-0052

10-0082