

I. APPLICATIONS

- 1. Cell centrifugation
- 2. Molecular biology
- 3. Sample storage
- 4. Centrifugation of precipitates

BIOLOGIX

II. FEATURES

- 1. Materials: Strong, medical-grade polypropylene, autoclavable
- 2. Chemically resistant to alcohols and mild organic solvents
- 3. Transparent wall permits easy viewing of tube contents
- 4. Working temperature: stable from -20°C to 121°C (Autoclavable)
- 5. Tested to withstand centrifugation of 20,000 RCF in a fully supported rotor
- 6. Certified RNase & DNase Free

III. TUBE DESIGN

- 1. Conical bottom design fits most standard floor model and table-top centrifuges
- 2. Molded graduations in 0.5ml increments
- 3. Clearly visible with either clear or turbid samples
- 4. Available in 4 colors: clear, red, yellow, and blue

IV. CAP DESIGN

1. Attached cap

V. PACKAGING

1. Bulk packaged, 500 tubes/pack and 10 packs/case

VI. DESCRIPTION

Item	Volumn	Tooth	Color	Packaging	Sterile
80-1500-S	1.5ml	No	Clear	Bulk, 500 tubes/pack, 10 packs/case	Yes
80-1500	1.5ml	No	Clear	Bulk, 500 tubes/pack, 10 packs/case	No
80-1501	1.5ml	No	Red	Bulk, 500 tubes/pack, 10 packs/case	No
80-1502	1.5ml	No	Yellow	Bulk, 500 tubes/pack, 10 packs/case	No
80-1503	1.5ml	No	Blue	Bulk, 500 tubes/pack, 10 packs/case	No
80-0500	0.5ml	Yes	Clear	Bulk, 500 tubes/pack, 10 packs/case	No
80-0015	1.5ml	Yes	Clear	Bulk, 500 tubes/pack, 10 packs/case	No
80-0020	2.0ml	Yes	Clear	Bulk, 500 tubes/pack, 10 packs/case	No

VII. STERILITY

1.N/A

Tel: +86-531-67802668 Fax: +86-531-67803768 www.BiologixGroup.com



VIII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

IX. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

X. PRODUCT PARAMETER

Dimensions=mm

80-1500 Tube 80-1500-S Tube


A = O.D. 13.0

B = O.D. 10.8

C = 39.3

D = O.D. 9.0

E = 1.5



80-1500/80-1500-S

80-0015 Tube


A = O.D. 13.0

B = O.D. 10.8

C = 39.3

D = O.D. 9.0

E = 1.5



Fax: +86-531-67803768

80-0015

Document Number: **BD-T0036**

Document Date: May 2014

Page 3/**4**

80-0500 Tube

Λ		\sim	. 1	\neg	\sim	_
Δ	=				9	5

B = O.D. 7.7


C = 30.4


D = O.D. 7.2

E = 4.4

F= 15.4

G = 31.7





80-0500

80-0020 Tube

A = O.D. 11.9

B = O.D. 9.0


C = 40.6


D = O.D. 8.9

E = 5.2

F= 18.7

G = 41.8





80-0020

Fax: +86-531-67803768



Document Number: **BD-T0036**

Document Date: May 2014

Page 4**/4**

Tel: +86-531-67802668

www.BiologixGroup.com

Fax: +86-531-67803768