Document Number: BD-T0169

Document Date: **Dec. 31**st, **2022**

Page 1/**3**

Laboratory Waste Fluid Collection System

I. APPLICATIONS

- 1. Professional design enables the safe disposal of chemicals
- 2.Anti-static, no danger caused by the static electricity
- 3. Superior sealing performance, super efficient filter
- 4. Accessories can be customized
- 5. Optional safety tray, lockable inside, leaking alarm

II. FEATURES

Funnel-type waste fluid collection device includes: one 20L waste fluid container with inner film, one inner cover, one outer cover, one funnel with a cover, one PTFE stop valve, one liquid lever alarm, one set of efficient filer connector, and one Inorganic/ organic chemical material filter.



Tel: +86-531-67802668 Fax: +86-531-67803768 www.BiologixGroup.com



Document Number: Document Date: Page BD-T0169 Dec. 31st, 2022 2/3

- Filmed inner wall
- Connection type: Can be connected with fluid instruments of various brands Funnel type: Equipped with the valve funnel
- High sealing performance
- Equipped with super efficient filer with connector
- Optional: Liquid level alarm

III. Product Specifications

Volume	Material	Dimensions (L*W*H) mm	Color	Funnel's Material	Outer Diameter of Funnel (mm)	Funnel' s Volume (ml)
20L	HDPE	272*272*400	Clear	PE-HD	150	600
20L	HDPE	400*270*300	Clear	PE-HD	150	600

- 1. Customized service. 10L, 20L and 60L containers, laboratory intelligent management system, carts, and safety trays are available upon request.
- 2. 5 types of filters are available, which are organic gas high-efficiency filter, inorganic gas high-efficiency filter, acid gas high-efficiency filter, ammonia and methylamine filter and multi-purpose high-efficiency filter, and high-efficiency filter.

IV. Product Specifications

ITEM	DESCRIPTION	PCS PER CASE
05-2721	Funnel-type waste fluid collection device-super efficient, One 20L waste fluid container with inner film, one inner cover, one outer cover, one funnel with a cover, one set of efficient filer connector, one Inorganic or organic chemical material filter	1
05-2791	Super efficient filter, Inorganic chemical material filter	1
05-2792	Super efficient filter, Organic chemical material filter	1
05-2793	Liquid lever alarm, Contact liquid lever alarm	1
05-2797	Liquid carrying cart (standard), Adjustable size; one cart and tray	1
05-2798	Liquid carrying cart (intelligent), Adjustable size; one cart and tray; equipped with intelligent module	1

Tel: +86-531-67802668 Fax: +86-531-67803768 www.BiologixGroup.com



Document Number: Document Date: Page BD-T0169 Dec. 31st, 2022 3/3

V. Regulatory Compliance Mark



VI. Warranty

The company will provide free maintenance if the quality problems are caused by non-human beings within 12 months from the date of delivery of the product.

Misused or unsatisfactory operation and the damage caused by unauthorized repair or modification shall not be covered by warranty. Repair costs outside the warranty period should be paid for.

CONTACT US

Biologix Group Limited Tel: +86-531-67802668

Web: www.BiologixGroup.com Email: Inquiry@BiologixGroup.com